

CENTRIFUGES

FIXED ANGLE ROTOR

Analogue

M 415 - M 615 - M 815 - M 4812M

Capacity	4 x 15 ml 6 x 15 ml 8 x 15 ml 12 x 15 ml
Rotor type	Vertical 35° fixed angle
Speed control	Stepless electronic
Max. speed	5000 rpm
Max. RCF	2515 x g (3075 x g M 4812M)
Timer range	1-15 min
Fuse	1 Amp
Power rating	100 W
Safety measure	Lid safety switch
Lid lock	Available
Rotor material	Polycarbonate
Tube socket (godet) material	ABS (rigid plastic)
Casing material	Sheet steel finished in stoved enamel paint + ABS
Power supply	AC 190-240 V / 50-60 Hz
Overall dimensions (w x h x d)	29 x 25 x 36 cm
Net weight	9 kg
Gross weight	10 kg



4 x 15 ml 6 x 15 ml 8 x 15 ml 12 x 15 ml

