

CENTRIFUGES

FIXED ANGLE ROTOR

Digital

M 450P

Capacity	4 x 50 ml
Rotor type	Vertical 35° fixed angle
Max. speed	5000 rpm
Max. RCF	3075 x g
Speed setting increment	10 rpm
Speed setting range	500 - 5000 rpm
Timer	Starts when set speed achieved
Timer range	1-99 min. In 1 minute step
Continuous running mode	Available
Running by entering RCF value	Available
Working short periods	With pulse button
Number of programs	10
Fuse	6 Amps
Power rating	350 W
Safety measure	Lid safety switch
Lid lock	Available
External speed reading	Available
Rotor material	Polycarbonate
Tube socket (godet) material	ABS (rigid plastic)
Casing material	Sheet steel finished in stoved enamel paint + ABS
Motor type	Maintenance free, brushless motor
Power supply	AC 190-240 V / 50-60 Hz
Overall dimensions (w x h x d)	33 x 30 x 45 cm
Net weight	16 kg
Gross weight	18 kg

4.3" TFT Touch Screen Display

Speed RCF and time control, acceleration / deceleration time

