

VORTEX MIXER

VORTEX MIXER

M 16

| | |
|--|--|
| Mixing method | <i>Vortex</i> |
| Mixing stroke height | <i>3 mm.</i> |
| Rotating speed range (cycle/min.) | <i>500 - 2500</i> |
| Speed control | <i>Non-gradual setting</i> |
| Power rating | <i>150 W</i> |
| Fuse | <i>1 Amp.</i> |
| Casing material | <i>Sheet steel finished in stoved enamel paint</i> |
| Power supply | <i>AC 190 - 240 V – 50/60 Hz</i> |
| Overall dimensions (w x h x d) | <i>18x 21 x 21 cm.</i> |
| Net weight | <i>7 kg.</i> |
| Gross weight | <i>7.5 kg.</i> |

