

SAND BATHS

SAND BATH

30 x 30 cm

M 3030K

Temperature working range	50 °C to 300 °C
Temperature control	Digital thermostat
Thermostat setting increments	in 1 °C steps
Thermostat working accuracy	± 1 °C
Power rating	1500 W
Fuse	16 Amps.
Casing material	Sheet steel finished in stoved enamel paint
Heating plate material	Normalized cast iron
Sand pool material	Stainless steel
Sand pool depth	10 cm
Heating plate dimensions	30 x 30 cm
Power supply	AC 190-240 V / 50-60 Hz
Overall dimensions (w x h x d)	30 x 24 x 30 cm
Net weight	15 kg
Gross weight	16 kg

