

CENTRIFUGES

FIXED ANGLE ROTOR

Digital

M 415 - M 615 - M 815 - M 4812P

Capacity	4 x 15 ml 6 x 15 ml 8 x 15 ml 12 x 15 ml.
Rotor type	Vertical 35° fixed angle
System control	PID digital microprocessor control
LCD Display	Displays speed RCF and time control
Max. speed	5000 rpm
Max. RCF	2515 x g (3075 x g M 4812P)
Speed setting increment	10 rpm
Speed setting range	500 - 5000 rpm
Timer	Starts when set speed achieved
Timer range	1-99 min. In 1 minute step
Continuous running mode	Available
Running by entering RCF value	Available
Process finalization signal	Audible signal available
Working short periods	With pulse button
Number of programs	10
Fuse	1 Amps
Power rating	100 W
Safety measure	Lid safety switch
Lid lock	Lid inter lock that prevent lid opening during centrifugation
Rotor material	Polycarbonate
Tube socket (godet) material	ABS (Rigid Plastic)
Casing material	Sheet steel finished in stoved enamel paint + ABS
Power supply	AC 190-240 V / 50-60 Hz
Overall dimensions (w x h x d)	29 x 25 x 36 cm
Net weight	9 kg
Packed weight	10 kg



4 x 15 ml 6 x 15 ml 8 x 15 ml 12 x 15 ml