

ELECTROSURGICAL UNITS

H- SERIES

400W Monopolar & Bi-polar

M 400

Cutting mode	<i>Available</i>
Mixed cutting-coagulation mode	<i>Available</i>
Coagulation mode	<i>Available</i>
RF maximum cutting power	<i>400 W</i>
RF maximum coagulation power	<i>100 W</i>
RF spray coagulation power	<i>80 W</i>
Bi-polar coagulation power	<i>80 W</i>
Bi-polar maximum cutting power	<i>100 W</i>
Power activation	<i>Button controlled pencil or foot switch</i>
Cutting coagulation impedance	<i>300Ω</i>
Spray coagulation impedance	<i>1000Ω</i>
Reusable electrode	<i>8 pcs</i>
Power rating	<i>400 W</i>
Fuse	<i>3 Amps</i>
Casing material	<i>Sheet steel finished in stoved enamel paint</i>
Power supply	<i>AC 220 V – 50/60 Hz</i>
Overall dimensions (w x h x d)	<i>36 x 16 x 38 cm</i>
Net weight	<i>10.5 kg</i>
Gross weight	<i>11.5 kg</i>

