

CENTRIFUGES

MICROHAEMATOCRIT CENTRIFUGE

M 19

Capacity (Capillary tubes)	24 pcs x (Ø2 x 75 mm)
Timer range	1 - 15 min
Max. speed	10000 rpm
Max. RCF	9000 x g
Safety measure	Lid safety switch
Lid lock	Available
Fuse	1 Amp
Power rating	100 W
Reading scale capacity	0-100 %
Reading scale measurement limits	45-70 mm
Tube holding disc material	Acetal resin
Tube holding disc cover material	Stainless steel
Casing material	Sheet steel finished in stoved enamel paint + ABS
Power supply	AC 190-240 V / 50-60 Hz
Overall dimensions (w x h x d)	29 x 25 x 36 cm
Net weight	9 kg
Gross weight	10 kg



24 pcs x (Ø2 x 75 mm)