

ELECTROSURGICAL UNITS

ELECTROSURGICAL UNITS

80W Monopolar & Bi-polar

M 20-80

| | |
|--------------------------------|---|
| Cutting mode | Available |
| Mixed cutting-coagulation mode | Available |
| Coagulation mode | Available |
| Bi-polar coagulation mode | Available |
| RF maximum cutting power | 80 W |
| RF maximum coagulation power | 40 W |
| RF spray coagulation power | 40 W |
| Bi-polar coagulation power | 16 W |
| Power activation | With foot switch |
| Cutting coagulation impedance | 300Ω |
| Spray coagulation impedance | 1000Ω |
| Reusable electrode | 8 pcs |
| Power rating | 120 W |
| Fuse | 2 Amps |
| Casing material | Sheet steel finished in stoved enamel paint |
| Power supply | AC 220 V – 50/60 Hz |
| Overall dimensions (w x h x d) | 30 x 12 x 24 cm |
| Net weight | 7.5 kg |
| Gross weight | 8 kg |

