

# WATER STILL

## WATER STILL

With Stainless Steel Chamber

### M 4

<b>Distillation water capacity</b>	4 lt / hour
<b>Cooling water consumption (max.)</b>	40 lt / hour
<b>Cooling water output temperature</b>	50 °C
<b>Heater type</b>	Stainless steel tube resistance
<b>Safety measure</b>	Water level control with thermostat and electrode
<b>Casing material</b>	Stainless steel
<b>Interior material</b>	Stainless steel
<b>Connection materials</b>	PVC and silicon hose
<b>Fuse (automatic)</b>	25 Amps.
<b>Power rating</b>	3000 W
<b>Monitoring influent pressure</b>	With analog manometer
<b>Power supply</b>	AC 190-240 V / 50-60 Hz
<b>Overall dimensions (w x h x d)</b>	65 x 52 x 18 cm.
<b>Net weight</b>	18 kg.
<b>Gross weight</b>	20 kg.



\*\*\*Optional Internal tank 7 lt