

AUTOCLAVES

AUTOCLAVE

Class - B

SM 18 B

Capacity	18 lt
Sterilizing Temperature	121 °C / 134 °C
Water Tank	Open-Upward Water Tank
Door Lock	Double-Controlled
Printer	Data Logger Printer
Dry Procedure	Vacuum Drying
Display	LED, English
Tests	Bowie - Dick & Vacuum Test
Safety Features	Pressure Release Valve and Computerized Diagnosing Including Alarm System
Capacity of Water Tank	Distilled water 3 lt / waste water 2,5 lt
Casing material	Sheet stainless steel finished in stoved enamel paint
Chamber material	Stainless steel
Power Supply	AC 190-240 V / 50-60 Hz
Power Rating	2000 W / 10A
Instrument Containers	3 trays
Overall Size (w x h x d)	
Net Weight	

