

CENTRIFUGES

FIXED ANGLE ROTOR

Analogue

M 415 - M 615 - M 815 E

Capacity	4 x 15 ml 6 x 15 ml 8 x 15 ml
Rotor type	Vertical 35° fixed angle
Speed control	Stepless electronic
Max. speed	5000 rpm
Max. RCF	2515 x g
Fuse	1 Amp.
Power rating	100 W
Safety measure	Lid safety switch
Rotor material	Polycarbonate
Tube socket (godet) material	ABS (rigid plastic)
Casing material	Sheet steel finished in stoved enamel paint + ABS
Power supply	AC 190-240 V / 50-60 Hz
Overall dimensions (w x h x d)	22 x 28 x 24 cm
Net weight	5.5kg
Gross weight	6kg



4 x 15 ml 6 x 15 ml 8 x 15 ml