

CENTRIFUGES

FIXED ANGLE ROTOR

Analogue

M 615 E

| | |
|---------------------------------------|---|
| Capacity | 6 x 15 ml |
| Rotor type | Vertical 35° fixed angle |
| Speed control | Stepless electronic |
| Max. speed | 5000 rpm |
| Max. RCF | 2515 x g |
| Fuse | 1 Amp. |
| Power rating | 100 W |
| Safety measure | Lid safety switch |
| Rotor material | Polycarbonate |
| Tube socket (godet) material | ABS (rigid plastic) |
| Casing material | Sheet steel finished in stoved enamel paint + ABS |
| Power supply | AC 190-240 V / 50-60 Hz |
| Overall dimensions (w x h x d) | 22 x 28 x 24 cm |
| Net weight | 5.5kg |
| Gross weight | 6kg |



6 x 15 ml