

# CENTRIFUGES

## FIXED ANGLE ROTOR

Analogue

### M 615 M

<b>Capacity</b>	6 x 15 ml
<b>Rotor type</b>	Vertical 35° fixed angle
<b>Speed control</b>	Stepless electronic
<b>Max. speed</b>	5000 rpm
<b>Max. RCF</b>	3057 x g (3075 x g M 4812M)
<b>Timer range</b>	1-15 min
<b>Fuse</b>	1 Amp
<b>Power rating</b>	100 W
<b>Safety measure</b>	Lid safety switch
<b>Lid lock</b>	Available
<b>Rotor material</b>	Polycarbonate
<b>Tube socket (godet) material</b>	ABS (rigid plastic)
<b>Casing material</b>	Sheet steel finished in stoved enamel paint + ABS
<b>Power supply</b>	AC 190-240 V / 50-60 Hz
<b>Overall dimensions (w x h x d)</b>	29 x 25 x 36 cm
<b>Net weight</b>	9 kg
<b>Gross weight</b>	10 kg



6 x 15 ml

