

CENTRIFUGES

FIXED ANGLE ROTOR

Digital

M 2415P

Capacity	24 x 15 ml
Rotor type	Vertical 35° fixed angle
Max. speed	5000 rpm
Max. RCF	3057 x g
Speed setting increment	10 rpm
Speed setting range	500 - 5000 rpm
Timer	Starts when set speed achieved
Timer range	1-99 min. In 1 minute step
Continuous running mode	Available
Running by entering RCF value	Available
Working short periods	With pulse button
Number of programs	10
Fuse	6 Amps.
Power rating	350 W
Safety measure	Lid safety switch
Lid lock	Available
External speed reading	Available
Rotor material	Acetal resin
Tube socket (godet) material	ABS (rigid plastic)
Casing material	Sheet steel finished in stoved enamel paint + ABS
Inner Chamber Material	Sheet aluminum finished in stoved enamel paint
Chamber material	Stainless steel
Motor type	Maintenance free, brushless motor
Power supply	AC 190-240 V / 50-60 Hz
Overall dimensions (w x h x d)	37 x 34 x 55 cm
Net weight	22 kg
Gross weight	24 kg

4.3" TFT Touch Screen Display

Speed RCF and time control, acceleration / deceleration time

