

HOT PLATES

HOT PLATE

Ø15 cm

M 15

Temperature working range	50 °C to 300 °C
Temperature control	Thermostat with liquid expansion
Thermostat setting increments	in 5 °C steps
Thermostat working accuracy	± 5 °C
Power rating	1000 W
Fuse	6 Amps
Casing material	Sheet steel finished in stoved enamel paint
Heating plate material	Normalized cast iron
Heating plate dimension	Ø15 cm
Power supply	AC 190 / 240 V – 50/60 Hz
Overall dimensions (w x h x d)	18 x 26 x 9 cm
Net weight	2 kg
Gross weight	2.5 kg

