

HOT PLATES

HOT PLATE

40 x 60 cm

M 4060

Temperature working range	50 °C to 300 °C
Temperature control	Digital thermostat
Thermostat setting increments	in 1 °C steps
Thermostat working accuracy	± 1 °C
Power rating	3000 W
Fuse	25 Amps.
Casing material	Sheet steel finished in stoved enamel paint
Heating plate material	Teflon coated sheet aluminum
Heating plate dimensions	40 x 60 cm
Power supply	AC 190 - 240 V / 50-60 Hz
Overall dimensions (w x h x d)	60 x 14 x 40 cm
Net weight	19 kg
Gross weight	24 kg

