

# WATER STILL

## WATER STILL

With Stainless Steel Chamber

### M 4

**Distillation water capacity**

4 lt/hour

**Cooling water consumption (max.)**

40 lt/hour

**Cooling water output temperature**

50 °C

**Heater type**

Stainless steel tube resistance

**Safety measure**

Water level control, magnetic switch, leakage current fuse

**Casing material**

Stainless steel

**Interior material**

Stainless steel

**Connection materials**

PVC and silicon hose

**Water filter**

External

**Internal tank capacity**

8 lt.

**Fuse (automatic)**

25 Amps.

**Power rating**

3000 W

**Monitoring influent pressure**

With analog manometer

**Power supply**

AC 190-240 V / 50-60 Hz

**Overall dimensions (w x h x d)**

65 x 52 x 25 cm.

**Net weight**

18 kg.

**Gross weight**

20 kg.

### M 8

8 lt/hour

80 lt/hour

50 °C

Stainless steel tube resistance

Water level control, magnetic switch, leakage current fuse

Stainless steel

Stainless steel

PVC and silicon hose

External

16 lt.

3x25 Amps.

3x2 kW

With analog manometer

AC 380 V / 50-60 Hz

72 x 66 x 34 cm.

32 kg.

35 kg.

### M 12

12 lt/hour

120 lt/hour

50 °C

Stainless steel tube resistance

Water level control, magnetic switch, leakage current fuse

Stainless steel

Stainless steel

PVC and silicon hose

External

24 lt.

3x25 Amps.

3x2,4 kW

With analog manometer

AC 380 V / 50-60 Hz

72 x 66 x 34 cm.

32 kg.

35 kg.



 elektro-mag®

Ikitelli O.S.B. Demirciler Sitesi, B7 Blok No:153 Basaksehir, 34490 - Istanbul / TURKEY

Tel: +90 212 549 55 25 Fax: +90 212 549 52 09

export@elektro-mag.com www.elektro-mag.com

CE