

NEONATAL EQUIPMENTS

INFANT WARMER

ECO

M 308

Working temperature range	32 °C – 39 °C
Type of thermoelement	Pt – 100
Temperature control	PID Digital control
Heather type	Ceramic
Thermostat setting increments	in 0.1 °C steps
Thermostat working accuracy	± 0.1 °C
Thermometer display increments	in 0.1 °C steps
Safety measure	Safety thermostat with liquid expansion
Body material	Sheet steel and profile finished in stoved enamel paint
Illumination from above	LED Panel
Acrylic railing	Available
Number of casters	4 (2 equipped with brakes)
Number of drawers	2
Fuse	6 Amps.
Power rating	1000W
Power supply	AC 190-240 V / 50-60 Hz
Overall dimensions (w x h x d)	80 x 172 x 55 cm.
Net weight	65 kg.
Gross weight	85 kg.

