

# NEONATAL EQUIPMENTS

## PHOTOTHERAPY EQUIPMENT ( 360° )

LED

### M 304

**Phototherapy light source**  
**Phototherapy lamp life**  
**Control board**  
**Avialable language selection**  
**Timer**  
**Alarm at the end of selected period**  
**Temperature regulating of thermostat**  
**Lamp life cycle counter**  
**Temperature range**  
**Skin Temperature Group**  
**Total operating range**  
**Operating intensity selection**  
**Power rating**  
**Fuse**  
**Main body material**  
**Main supporting material**  
**Number of casters**  
**Power supply**  
**Netting mattress dimensions**  
**Overall dimensions (w x max.h x d)**  
**Total net weight**  
**Total gross weight**

*LED lamp*  
*50.000 hours*  
*4.3 " colored TFT touch screen display*  
*Turkish, English, French, Arabic, Russian*  
*Digital*  
*Available*  
*Available*  
*Available*  
*32-38 C*  
*Available*  
*0-12 hours*  
*%25, %50, %75, %100*  
*150 W*  
*2 Amp*  
*Sheet steel finished in stoved enamel paint*  
*Steel profile finished in stoved enamel paint*  
*4 ( 2 equipped with brakes )*  
*AC 190-240 V / 50-60 Hz*  
*40 x 65 cm*  
*56 x 145 x 75 cm*  
*35 kg*  
*45 kg*




 **elektro·mag®**

LABORATUVAR ALETLERİ SAN. VE TİC. A.Ş.

 1984

İkitelli O.S.B. Demirciler Sanayi Sitesi B7 Blok No:153 34490 - Basaksehir İSTANBUL

 0 212 5349 55 25 pbx

 0 212 549 52 09

 [export@elektro-mag.com](mailto:export@elektro-mag.com)

 [www.elektro-mag.com](http://www.elektro-mag.com)