


# NEONATAL EQUIPMENTS

## PHOTOTHERAPY EQUIPMENT

LED

### M 304

<b>Phototherapy light source</b>	LED lamp
<b>Phototherapy lamp life</b>	50.000 hours
<b>Control board</b>	3.5" colored TFT touch screen display
<b>Avialable language selection</b>	Turkish, English, French, Arabic, Russian
<b>Timer</b>	Digital
<b>Alarm at the end of selected period</b>	Available (audible)
<b>Lamp life cycle</b>	Available
<b>Total operating range</b>	0-12 hours
<b>Operating intensity selection</b>	50, %75, %100
<b>Light power</b>	60 MW
<b>Light wavelength</b>	460 nm
<b>Power rating</b>	25 W
<b>Fuse</b>	2 Amp
<b>Main body material</b>	ABS plastic
<b>Main supporting material</b>	Steel profile finished in stoved enamel paint
<b>Number of casters</b>	4 (2 equipped with brakes)
<b>Power supply</b>	AC 190-240 V / 50-60 Hz
<b>Dimensions of phototherapy unit (w x h x d)</b>	25 x 8 x 45 cm
<b>Overall dimensions (w x max.h x d)</b>	50 x 190 x 80 cm
<b>Total net weight</b>	14 kg
<b>Total gross weight</b>	20 kg



 elektro·mag®

LABORATUVAR ALETLERİ SAN. VE TİC. A.Ş.

İkitelli O.S.B. Demirciler Sanayi Sitesi B7 Blok No:153 34490 - Başakşehir İSTANBUL

☎ 0 212 549 55 25 pbx



0 212 549 5209



satis@elektro-mag.com



www.elektro-mag.com

CE 1984

Certificate Number : 1984-MDD-18-509