

# ELECTROSURGICAL UNITS

## H- SERIES

400W Monopolar & Bi-polar

Cutting mode	Available
Mixed cutting-coagulation mode	Available
Coagulation mode	Available
Maximum cutting power	400 W
Maximum coagulation power	100 W
Spray coagulation power	70 W
Bi-polar coagulation power	70 W
Bi-polar maximum cutting power	100 W
Power activation	Button controlled pencil or foot switch
Cutting coagulation impedance	300Ω
Spray coagulation impedance	300Ω
Reusable electrode	6 pcs
Power rating	400 W
Fuse	3 Amps
Casing material	Sheet steel finished in stoved enamel paint+plastic
Power supply	AC 190 - 240 V / 50 - 60 Hz
Overall dimensions (w x d x h)	36 x 42 x 17 cm
Net weight	10.5 kg
Gross weight	11.5 kg



elektro·mag®

[www.elektro-mag.com](http://www.elektro-mag.com)

CE

1984

Certificate Number : 1984-MDD-18-509

MED