

# NEGATOSCOPES

## PEGA IV

<b>Dimensions (h x w)</b>	520 x 1620 mm
<b>Screen Size (h x w)</b>	430 x 1525 mm
<b>Power Consumption</b>	140 W
<b>Expected Lamp Life Time (Aprox.)</b>	100.000 Hours
<b>Frame</b>	Anodized Aluminum
<b>Net Weight</b>	12,3 kg



World's thinnest (2.4cm),  
most homogeneous LED negatoscope!  
Very low energy consumption  
Expected Lamp life time approx. 100.000 hours



LABORATUVAR ALETLERİ SAN. VE TİC. A.Ş.

