

# ELECTROSURGICAL UNITS

## RF ELECTROSURGICAL UNITS

120W Monopolar & Bi-polar

### M 20-80

Monopolar cut	120 W
Monopolar cut / coagulation	80 W
Monopolar Hemo	60 W
Monopolar Fulgurate	40 W
Bipolar	120 W
Monopolar cut impedance	500 $\Omega$
Monopolar cut / coagulation impedance	500 $\Omega$
Monopolar hemo impedance	500 $\Omega$
Monopolar fulgurate impedance	500 $\Omega$
Bipolar impedance	200 $\Omega$
Power activation	Button controlled pencil or foot switch
Hand manipulator electrode	12 pcs
Power rating	120 W
Fuse	2 Amps
Casing material	Sheet steel finished in stoved enamel paint+plastic
Power supply	AC 190-240 V / 50-60 Hz
Overall dimensions (w x h x d)	30 x 12 x 24 cm
Net weight	7.5 kg
Gross weight	8 kg

