

# NEONATAL EQUIPMENTS

## PHOTOTHERAPY EQUIPMENT

LED

### M 304

**Phototherapy light source**

**Phototherapy lamp life**

**Control board**

**Available language selection**

**Timer**

**Alarm at the end of selected period**

**Lamp life cycle**

**Total operating range**

**Operating intensity selection**

**Light power**

**Light wavelength**

**Maximum light intensity (from 40 ± 1 cm distance)**

**Effective area (effective surface area)**

**Operation mode**

**Power rating**

**Fuse**

**Main body material**

**Main supporting material**

**Number of casters**

**Power supply**

**Dimensions of phototherapy unit (w x h x d)**

**Overall dimensions (w x max.h x d)**

**Total net weight**

**Total gross weight**

Blue fluorescent LED lamp

70.000 hours

3.5" colored TFT touch screen display

Turkish, English, French, Arabic, Russian

Digital

Available (audible)

Available

0-12 hours

4 Mode (%25, %50, %75, %100)

60 MW

460 nm

60  $\mu\text{w}/\text{cm}^2/\text{nm}1$  ( $\pm 10$ )

40 cm x 20 cm (120  $\mu\text{w}/\text{cm}^2/\text{nm}1$  ( $\pm 5$  cm)

3 mode (low, medium, high)

25 W

2 Amp

ABS plastic

Steel profile finished in stoved enamel paint

4 (2 equipped with brakes)

AC 190-240 V / 50-60 Hz

25 x 8 x 45 cm

50 x 190 x 80 cm

14 kg

20 kg

20 kg



 elektro·mag®

LABORATUVAR ALETLERİ SAN. VE TİC. A.Ş.

İkitelli O.S.B. Demirciler Sanayi Sitesi B7 Blok No:153 34490 - Başakşehir İSTANBUL

☎ 0 212 549 55 25 pbx



0 212 549 5209



satis@elektro-mag.com



www.elektro-mag.com

CE 1984

Certificate Number : 1984-MDD-18-509