

NEONATAL EQUIPMENTS

PHOTOTHERAPY EQUIPMENT (360°)

LED

M 304

Phototherapy light source	LED lamp
Phototherapy lamp life	50.000 hours
Control board	7 " colored TFT touch screen display
Avialable language selection	Turkish, English, French, Arabic, Russian
Timer	Digital
Alarm at the end of selected period	Available
Temperature regulating of thermostat	Available
Lamp life cycle counter	Available
Temperature range	32-38 C
Skin Temperature Group	Available
Total operating range	0-12 hours
Operating intensity selection	%25, %50, %75, %100
Maximum light intensity (from 40 ±1cm distance)	120 µw/cm ² /nm1 (±%10)
Minimum light intensity (from 40 ±1cm distance)	30 µw/cm ² /nm1 (±%10)
Power rating	150 W
Fuse	2 Amp
Main body material	Sheet steel finished in stoved enamel paint
Main supporting material	Steel profile finished in stoved enamel paint
Number of casters	4 (2 equipped with brakes)
Power supply	AC 190-240 V / 50-60 Hz
Netting mattress dimensions	40 x 65 cm
Overall dimensions (w x max.h x d)	56 x 145 x 75 cm
Total net weight	35 kg
Total gross weight	45 kg



 elektro·mag®

LABORATUVAR ALETLERİ SAN. VE TİC. A.Ş.

CE 1984

İkitelli O.S.B. Demirciler Sanayi Sitesi B7 Blok No:153 34490 - Basaksehir İSTANBUL

☎ 0 212 5349 55 25 pbx

☎ 0 212 549 52 09

✉ export@elektro-mag.com

🌐 www.elektro-mag.com