

FUME HOODS

FUME HOOD

With Cabinet

EM 901

Fan suction capacity	1750 m ³ /hour
Motor speed control	Available
Sliding cabin glass	Laminated safety glass 6 mm
Gas and water valves	Available
Gas and water nozzles	Available
Power outlet	2 pcs
Cabin illuminating	Available
Stainless steel sink	Available
Cabinet with casters	Available
Working surface material	Acid proof laminate
Casing material	Galvanized steel finished in stoved enamel paint
Fuse	2 x 25 Amps.
Power rating	100 W
Power supply	AC 190-240 V / 50-60 Hz
Working cabin dimensions (w x h x d)	88 x 90 x 75 cm.
Overall dimensions (w x h x d)	90 x 225 x 88 cm.
Net weight	120 kg.
Gross weight	160 kg.

