

# FUME HOODS

## FUME HOOD

Bench top

### EM 902

<b>Fan suction capacity</b>	1750 m <sup>3</sup> /hour
<b>Motor speed control</b>	Available
<b>Sliding cabin glass</b>	Laminated safety glass 6 mm
<b>Gas and water valves</b>	Available
<b>Gas and water nozzles</b>	Available
<b>Power outlet</b>	2 pcs
<b>Cabin illuminating</b>	Available
<b>Stainless steel sink</b>	Available
<b>Working surface material</b>	Acid proof laminate
<b>Casing material</b>	Galvanized steel finished in stoved enamel paint
<b>Fuse</b>	2 x 25 Amps.
<b>Power rating</b>	100 W
<b>Power supply</b>	AC 190-240 V / 50-60 Hz
<b>Working cabin dimensions (w x h x d)</b>	90 x 118 x 75 cm.
<b>Overall dimensions (w x h x d)</b>	120 x 152 x 88 cm.
<b>Net weight</b>	65 kg.
<b>Gross weight</b>	100 kg.

