

# NEONATAL EQUIPMENTS

## PHOTOTHERAPY EQUIPMENT ( 360° )

LED

### M 304

<b>Phototherapy light source</b>	Blue fluorescent LED lamp
<b>Phototherapy lamp life</b>	70.000 hours
<b>Control board</b>	7 " colored TFT touch screen display
<b>Avialable language selection</b>	Turkish, English, French, Arabic
<b>Timer</b>	Digital
<b>Light wavelength</b>	460 nm
<b>Alarm at the end of selected period</b>	Available
<b>Lower / upper skin temperature threshold</b>	Available
<b>Lower / upper air temperature threshold</b>	Available
<b>Temperature regulating of thermostat</b>	Available
<b>Lamp life cycle counter</b>	Available
<b>Baby mattress</b>	Sliding out of the tunnel ( with alarm )
<b>Temperature range</b>	32-38 °C
<b>Separately or together operable lower and upper light source sections</b>	Available
<b>Skin Temperature Group</b>	Available
<b>Total operating range</b>	0-12 hours
<b>Operating intensity selection</b>	%25, %50, %75, %100
<b>Maximum light intensity (from 40 ±1cm distance)</b>	120 μw/cm <sup>2</sup> /nm1 (±%10)
<b>Minimum light intensity (from 40 ±1cm distance)</b>	30 μw/cm <sup>2</sup> /nm1 (±%10)
<b>Power rating</b>	150 W
<b>Fuse</b>	2 Amp
<b>Main body material</b>	Sheet steel finished in stoved enamel paint
<b>Main supporting material</b>	Steel profile finished in stoved enamel paint
<b>Number of casters</b>	4 (2 equipped with brakes)
<b>Power supply</b>	AC 190-240 V / 50-60 Hz
<b>Netting mattress dimensions</b>	40 x 65 cm
<b>Overall dimensions (w x max.h x d )</b>	56 x 145 x 75 cm
<b>Total net weight</b>	35 kg
<b>Total gross weight</b>	45 kg



 elektro·mag®

LABORATUVAR ALETLERİ SAN. VE TİC. A.Ş.

CE 1984

İkitelli O.S.B. Demirciler Sanayi Sitesi B7 Blok No:153 34490 - Basaksehir İSTANBUL

☎ 0 212 5349 55 25 pbx

☎ 0 212 549 52 09

✉ [export@elektro-mag.com](mailto:export@elektro-mag.com)

🌐 [www.elektro-mag.com](http://www.elektro-mag.com)