

WATER BATHS

WATER BATH

With Stainless Steel Chamber

M 920KP

Temperature working range	Ambient Temperature + 5 °C to 100 °C
Capacity (Tank volume)	30 lt
Type of thermoelement	NTC
Temperature control	PID Digital control
Load switching element	SSR (Solid State Relay)
Thermostat setting increments	in 0.1 °C steps
Thermostat working accuracy	± 0.1 °C
Thermometer display increments	in 0.1 °C steps
Tank temperature accuracy	± 0.5 °C
Postpone / Process time	1 min – 99 hours 59 min
Timing begins	When reached set temperature value
Mixing	Available
Fuse	6 Amps
Tank material	Stainless steel
Casing material	Stainless steel
Tube capacity	96 pcs
Number of tube support	4 pcs
Power rating	1250W
Power supply	AC 190-240 V / 50-60 Hz
Chamber dimensions (w x d x h)	50 x 30 x 20 cm
Overall dimensions (w x d x h)	63 x 35 x 26 cm
Net weight	18 kg
Gross weight	20 kg

