

AUTOCLAVES

SM 23B

CLASS - B

TEKNİK ÖZELLİKLER

Fuctional state	<i>New</i>
Sterilization method	<i>Steam temperature and pressure sterilization</i>
Stationary equipment	<i>Floor standing , front loading , single door</i>
Sterilization material	<i>Unwrapped instrument, packaged tool, plastic products, fibrous material</i>
Sterilization procedure	<i>Sterilization preparation by 3-folsd alternation of positive and negative pressure and sterilization by selection of temperature and pressure</i>
Maximum temperature	<i>139 °C</i>
Maximum pressure	<i>0.25 Mpa / 2.5 bar</i>
Minimum pressure	<i>-0.9 bar</i>
Minimum sterilization volume	<i>23 lt</i>
Shelves for tools and materials	<i>At least 3 pieces</i>
Protection against excess temperature and pressure	<i>Automatic protection: over temperature/pressure protection : low water / anti dry burning ; safety door lock</i>
Chamber volume	<i>Ø250 x 470</i>
Overall dimensions (w x h x d)	<i>505 x 710 x 450</i>
Packing dimensions (w x h x d)	<i>760 x 560 x 510</i>
Voltage	<i>110 V / 220 v , 60 Hz / 50 Hz</i>
Power	<i>2000 W</i>
Net weight	<i>58 kg</i>
Gross weight	<i>65 kg</i>



 elektro-mag®

LABORATUVAR ALETLERİ SAN. VE TİC. A.Ş.

Turgut Özal Caddesi Karagül İş Merkezi No:84 / 5 34093 - Fındıkzade İSTANBUL

☎ 0 212 534 8026-27

☎ 0 212 635 4060

✉ satis@elektro-mag.com

🌐 www.elektro-mag.com