

AUTOCLAVES

AUTOCLAVE

Class - B

SM 50 B

Fuctional state	<i>New</i>
Sterilization method	<i>Steam temperature and pressure sterilization</i>
Stationary equipment	<i>Floor standing, front loading, single door</i>
Sterilization material	<i>Unwrapped instrument, packaged tool, plastic products, fibrous material</i>
Sterilization procedure	<i>Sterilization preparation by 3-folsd alternation of positive and negative pressure and sterilization by selection of temperature and pressure</i>
Maximum temperature	<i>139 °C</i>
Maximum pressure	<i>0.25 Mpa / 2.5 bar</i>
Minimum pressure	<i>-0.9 bar</i>
Minimum sterilization volume	<i>50 lt</i>
Shelves for tools and materials	<i>At least 3 pieces</i>
Protection against excess temperature and pressure	<i>Automatic protection: over temperature/pressure protection : low water / anti dry burning ; safety door lock</i>
Chamber volume	<i>Ø316 x 640</i>
Overall dimensions (w x h x d)	<i>657 x 870 x 575</i>
Packing dimensions (w x h x d)	<i>980 x 750 x 655</i>
Voltage	<i>110 V / 220 v, 60 Hz / 50 Hz</i>
Power	<i>2800 W</i>
Net weight	<i>108 kg</i>
Gross weight	<i>142 kg</i>

