

SAFETY CABINETS

CLASS-II TYPE A / B

Operator, Product and Environmental Protections

Fan Capacity	2200 m ³ / hour
Motor Sound Level	< 60 db
Sliding Cabin Glass	Tempered 6 mm
Min. Inflow Air Speed	0,35 - 0,51 m/sn
Motor Speed Button	4.3" TFT Touch screen
Motor Control	4.3" TFT Touch screen
Electrical Outlet	2 pcs
U.V Light	1 pc (Uvc Type 240 - 260 nm)
Day Light	1 pc
Hepa Filter	2 pcs (99,999% Efficiency for Particles ≥0,3µm H14 class according to EN 1822)
Hepa Filter Life Timer	4.3" TFT Touch screen
Working Surface Material	Stainless Steel (AISI 304)
Casing Material	Galvanized Steel Finished in Electrostatic Powder Coating
Power Supply	AC 190-240 V / 50-60 Hz

Working Principles

Type A 1	%70 Exhaust , %30 Recirculate
Type A 2	%30 Exhaust , %70 Recirculate
Type B 1	%70 Exhaust and duct , %30 Recirculate
Type B 2	%100 Exhaust and duct



LABORATUVAR ALETLERİ SAN. VE TİC. A.Ş.